



PCR or LAMP

Food Safety Considerations when Choosing Molecular Detection Methods

Dr. Vikrant Dutta and Joy Dell'Aringa

Food microbiology pathogen detection technology is constantly evolving and improving for fast, efficient, and accurate analysis. Thanks to the wide commercialization of easy to use diagnostic kits, the end-user no longer needs a deep understanding of the intricacies of diagnostic chemistries to perform the analysis. However, when navigating the selection process in search of the technology that is best fit-for-purpose, it is critical to understand the key differences in principle of detection and how they can impact both operations and risk. Here, we will explore the difference between two broad categories of molecular pathogen detection: PCR and isothermal technologies such as LAMP.

PCR & LAMP Detection Chemistries: An Overview

PCR detection chemistries have come a long way from non-specific DNA-binding dyes like SYBR Green, to highly precise sequence-specific molecular probes. The efficiency of the Real-Time PCR reaction today allows for the use of a variety of detection probes, the most popular being Dual-Labeled Fluorescent Probes such as FRET, TaqMan probes, and Molecular Beacon probes¹. The precision of these probes is showcased in their ability to distinguish allelic single-nucleotide polymorphisms (SNPs)^{2,3}. The most prevalent isothermal chemistry: Loop-Mediated Isothermal Amplification (LAMP), typically does not use molecular probes due to the lack of structure and formation consistency in its amplified products. As a result, LAMP mostly relies on detection through non-specific signal generation like ATP bioluminescence or non-specific dyes. In theory, this could come from specific and non-specific amplification events. This also makes LAMP inept to detect the allelic polymorphisms, which in some cases are critical to detecting crucial variations, like between close species, and within serotypes. In the end, the detection chemistries are only as good as the amplified products.

Key Take Aways:

- PCR technology has improved greatly in detection efficiencies via target specific probes
- LAMP technology typically does not utilize specific molecular probes, but instead relies on indirect signal generation
- Target specific probes ensures signal from specific amplification events only
- Indirect signal can come from specific and non-specific amplification events, which can lead to a reduced specificity and inability to detect in certain cases

PCR & LAMP: Amplification Strategies

Food safety pathogen detection protocols aim to find the single cell of target organism lurking in a relatively large sample. In order to achieve detection, molecular technologies utilize amplification strategies to increase the concentration of target DNA to a detectable level. Nucleic acid amplifications in both PCR and isothermal technologies start by making a variety of amplified products. These products include



PCR or LAMP: Food Safety Considerations when Choosing Molecular Detection Methods

Non-Specific Amplifications (NSA), and specific (target) amplifications^{4,5,6,7}. Ideally, the concentration of the desired target amplified product increases over time to levels above NSA where the detection chemistries are able to provide a detectable signal from the desired amplified product (target). Various reaction components such as: target DNA concentration, polymerase, buffers, and primers play a defining role in maintaining the progressive amplification dynamics, and thereby act as core contributors to the robustness of the reaction. However, none play a more crucial contribution to the success of a reaction than temperature. Herein lies a key difference between the fundamentals of PCR and Isothermal amplification technologies.

Key Take-Aways:

- PCR and LAMP both make a variety of amplification products – Non-Specific (NSA) and Specific (target)
- Ideally, target products increase above the levels of NSA to reach a reliable detectable signal
- A variety of factors contribute to the overall robustness of the reaction

What is the difference between PCR and Isothermal Detection Technologies?

A key foundational difference between the two technologies lies in the utilization of the thermal profiles. PCR utilizes thermocycling, while isothermal does not. This difference is the tether around how the different amplification chemistries work. In PCR, the cyclical denaturation of DNA during thermocycling separates all dimers (specific and non-specific). As the reaction progresses, this leads to frequent correction of the amplification dynamics away from the NSA and favors amplification of the desired target amplifications. Isothermal chemistries do not have the ability to correct the NSA through thermocycling, so it must rely on alternate mechanisms to achieve the same result. For example, LAMP utilizes “nested” primers where the

primer sequences outside the target region are used to create early amplification products. These are subsequently used as a template for the desired target amplifications. The presence of these extra primers, along with the diverse amplified structures formed during the LAMP reaction, creates many more opportunities for NSA production^{8,9,5}. This causes a less controlled and inefficient amplification, and is perhaps why the preheating of the DNA prior to the LAMP has shown to increase the LAMP sensitivity^{10,11}. To the end user, this inefficiency can manifest itself in various ways such as restricted multiplexing, lack of internal amplification control, complex assay design, tedious sample prep methods, and increased chance for inaccurate results (i.e. false positives and false negatives)¹². Scientific literature does provide a fair amount of evidence that, under controlled conditions, the isothermal amplification reaction can provide equivalent results to PCR. Isothermal chemistries also usually require simplified instruments and thereby can present interesting opportunities in non-conventional test environments with simple and predictable matrices. This likely explains the early footing of isothermal technologies in the clinical test environment as a “point of care test” (POCT) alternative. However, it must also be noted that recently PCR has also been adapted and successfully commercialized for the POCT format^{13,14}.

Key Take-Aways:

- PCR utilizes thermocycling, Isothermal does not
- In PCR, thermocycling allows for the reaction to favor the target amplification over the NSA
- LAMP must rely on alternate mechanisms to correct for NSA and these mechanisms lead to a less controlled and therefore inefficient amplification
- Under controlled conditions, isothermal technology can provide equivalent results to PCR
- Low instrumentation requirements make Isothermal technologies interesting for non-conventional test environments (i.e. POCT); however, PCR has also been recently adapted as a POCT



PCR or LAMP: Food Safety Considerations when Choosing Molecular Detection Methods

Internal Amplification Controls in Molecular Pathogen Detection Technologies: The Value & The Challenges

The purpose of an Internal Amplification Control (IAC) is to provide an indication of the efficacy of the test reaction chemistry. The closer the IAC is to the target DNA sequence, the better view into the inner workings of each reaction. For food microbiology testing, the role of the Internal amplification control is more important now than ever before. Driven by regulations, industry self-accountability, and brand protection initiatives, more and more food laboratories are testing diverse product types with novel and innovated formulations and ingredients. IAC capability not only helps with troubleshooting, but it also allows for a more confident adoption of the technology for new and diverse food and environmental matrices.

Over the years, PCR has progressively developed into a robust and efficient technology that can provide a dynamic IAC, giving the end user a direct look into the compatibility of the test matrix within the PCR reaction. From a single reaction, we can now make a qualitative assessment of whether the crude DNA prep from a matrix undergoing testing is working with this PCR or if it is inhibiting the reaction. With legacy technologies, including the older generation PCR's, we were limited to an "it-did-not-work" scenario, leaving the end user blind to any insights into the reason. Since isothermal chemistries typically do not have an IAC, the end user is vulnerable to false results. Even when isothermal chemistries such as Nicking Enzyme Amplification Reaction -(NEAR) can provide IAC, they typically do not mimic the target reaction and, therefore, are not a direct indicator of the reaction dynamics. This limits the end user back to the "it-did-not-work" scenario. LAMP technology attempts to mitigate the absence of IAC by performing a separate, external, reaction with each test matrix. This strategy leaves the final result vulnerable to a number of factors that are otherwise non-existent for IAC: sampling variations, reagent and machine anomalies, and user error. External control approaches also have a notable impact to the end user, as the burden to demonstrate

fit-for-purpose of the method for even the smallest matrix composition change increases both validation and verification activities, which can have a notable financial impact to the laboratory.

There are a few reasons why IAC incorporation is not always plausible for isothermal technologies such as LAMP. First, inefficient, less-controlled amplification reactions leave little room for reliable and meaningful supplementary reactions, like the ones required for IAC. Second, the lack of consistent amplified products make it much more difficult to pinpoint a DNA structure that can be dependably used as an IAC. Third, lack of specific detection mechanisms makes it hard to distinguish signal from the target versus the IAC reaction.

Key Take-Aways:

- Internal amplification controls (IAC) are critical for the food industry due to complex and ever-changing matrix formulations
- IAC is useful for troubleshooting, optimizing assay performance, and adapt test for novel matrices
- PCR has evolved to provide dynamic IAC, leading to increased confidence in results
- LAMP is not able to utilize IAC due to the nature of the amplification products, reaction efficiency, and lack of specific detection mechanisms

How Does Sample Preparation Differ Between PCR & LAMP?

Post-incubation sample preparation for both PCR and Isothermal chemistries must meet two criteria: make the DNA available for amplification and remove any components that may interfere with the chemistry of the downstream detection method. In the context of food pathogen testing, the extraction of DNA is rarely an issue; it's the removal of the interfering PCR inhibitors that typically pose the greatest challenge. Much work has been



PCR or LAMP: Food Safety Considerations when Choosing Molecular Detection Methods

done to document the nature of PCR inhibitors and their mechanism of action¹⁵. This information has led to a co-evolution of PCR reagents and assay design parameters over the years. Today, we have PCR chemistries that can work with essentially “universal” sample prep method^{16,17,18,19,20,21}. This has done nothing short of revolutionize the workflow while ensuring optimal assay performance with minimal sample manipulations. Isothermal chemistries, including LAMP, continue to lag behind in this area. The lack of controlled reaction makes it aberrant to work with various kinds of food matrices such that they may work robustly with one matrix and would not with the next. More often than desired, sample prep methods have to be adjusted for compatibility, or contingencies are built into the “unified” sample prep to account for these aberrations, making them more complex than desired. For the time being, the isothermal sample prep workflow is expected to remain multi-step and nuanced. A natural consequence of this for the end user translates to inconsistent results.

Key Take-Aways:

- Removal of interfering components is essential in sample preparation
- PCR has evolved dramatically to address inhibitors and now offers robust “universal” solutions for easy and efficient sample preparation
- Isothermal is still developing in this area, and currently requires many adjustments to optimize sample preparation based on the matrices tested

Food Industry Applications - Which is Best Suited?

Among the molecular methods, both PCR and Isothermal technologies have a place and will continue to shape the food safety industry. The evolutionary genesis of these technologies can provide clues into the specific suitability for the different stakeholders within the food industry. PCR has come a long way. Over the years, many users

have reported on its strengths and shortcomings and a continuous development path has led to reliable technology that continues to be a benchmark for the molecular detection technologies. Other adaptable and advanced capabilities, like dynamic IAC inclusion, multiplexing, proven compatibility to peripherals like robotics for high throughput testing, adaptation for POCT, and easy access will continue to make PCR a favorable solution in the foreseeable future. Isothermal amplification technologies, including LAMP, provide interesting opportunities in non-conventional test environments, presenting viable alternatives for non-complex matrices, especially in the low-resource environments. However, with a deeper look, it is not hard to see that the isothermal chemistries are still grappling with the similar technical issues that beleaguered PCR more than a decade ago. From a tangential view, with the recent advancements and shift-in-focus to the “-omics” technologies, it can be argued that the need to develop or improve on other PCR-like technologies is likely going to take a back seat.

Key Take-Aways:

- PCR has the advantage of decades of evolution and improvement allowing for IAC, multiplexing, and ease of use – which for the time being, favor the risk-based industry
- LAMP has the opportunity to mimic the advances of PCR in the coming years

As new molecular technologies are ever emerging, the food safety industry is taking a close look at how the applications will impact the safety and integrity of our food supply. Standard consideration factors such as time to result, validations, sensitivity, and specificity will continue to lead the discussion. However, a deeper look at the chemistry bears scrutiny as well. Amplification efficiency, direct vs indirect signal detection, and internal amplification control inclusion are also being examined by industry leaders as molecular pathogen technologies continue to evolve.



PCR or LAMP: Food Safety Considerations when Choosing Molecular Detection Methods

Peer Review Acknowledgement

Technical review provided by*:

- Benjamin Bastin, Q Laboratories
- Benjamin Howard, Certified Laboratories
- Dr. Leslie Thompson-Strehlow, SGS Vanguard Sciences, Inc.

References

1. Navarro et al. (2015). Real-time PCR detection chemistry. *Clin Chim Acta.*, 439, 231-50.
2. Malkki Petersdorf. (2012). Genotyping of single nucleotide polymorphisms by 5' nuclease allelic discrimination. *Methods Mol Biol.*, 882, 173-82.
3. Sobrino et al. (2005). SNPs in forensic genetics: review on SNP typing methodologies. *Forensic Sci Int.*, 154(2-3), 181-94.
4. Ruiz-Villalba, A., et al. (2017). Amplification of nonspecific products in quantitative polymerase chain reactions (qPCR). *Biomol Detect Quantif.*, 14, 7-18.
5. Senarathet al. (2014). Non Specific Amplification with the LAMP Technique in the Diagnosis of Tuberculosis in Sri Lankan Settings. *JTR*, 2, 168-172.
6. Wang. (2016). Effect of internal primer-template mismatches on loop mediated isothermal amplification. *Biotechnology & Biotechnological Equipment*, 30(2), 314-318.
7. Gadkar et al. (2018). Real-Time Detection and Monitoring of Loop Mediated Amplification (LAMP) Reaction Using Self-quenching and De-quenching Fluorogenic Probes. *Sci Rep*, 8(1), 5548
8. Kimura et al. (2011). Optimization of turn-back primers in isothermal amplification. *Nucleic Acids Res.*, 39, e59.
9. Lee et al. (2009). Detection of genetically modified organisms (GMOs) using isothermal amplification of target DNA sequences. *BMC Biotechnol*, 9, 7.
10. Ihira et al. (2008). Loop-mediated isothermal amplification for discriminating between human herpesvirus 6A and B. *J Virol Methods.*, 154, 223-225.
11. Tomita et al. (2008). Loop-mediated isothermal amplification (LAMP) of gene sequences and simple visual detection of products. *Nat Protoc*, 877-882.
12. Wong, Yien-Ping & Othman, Shuhaidah & Lau, Yee-Ling & Son, Radu & Chee, H-Y. (2017). Loop Mediated Isothermal Amplification (LAMP): A Versatile Technique for Detection of Microorganisms. *Journal of Applied Microbiology*. 124. 10.1111/jam.13647.
13. Piralla et al. (2017). FilmArray(TM)GI panel performance for the diagnosis of acute gastroenteritis or hemorrhagic diarrhea. *BMC Microbiol*, 17(1), 111.
14. Buss et al. (2015). Multicenter evaluation of the BioFire FilmArray gastrointestinal panel for etiologic diagnosis of infectious gastroenteritis. *J Clin Microbiol*, 53(3), 915-25.
15. Schrader et al. (2012). PCR inhibitors - occurrence, properties and removal. *J. Appl Microbiol.*, 113(5), 1014-26.
16. Dutta et al. Matrix Extension Validation Study of the GENE-UP® Salmonella Assay for AOAC PTM. AOAC annual meeting, 2017, Atlanta, GA
17. Bird et al. An Independent Evaluation of the GENE-UP® EHEC Detection Method for the Detection of Non-O157 Shiga-toxin Producing *Escherichia coli* (STEC) and *Escherichia coli* O157:H7 in Foods. 2017, IAFF, Tampa, FL
18. Bird et al. An Independent Evaluation of the GENE-UP *Listeria monocytogenes* Assay for the Detection of *Listeria monocytogenes* in Foods. 2017, IAFF, Tampa, FL
19. Mathia et al. EN ISO 16140-2 Validation Study of the GENE-UP® PCR Method for the Detection of *Listeria* sp. in a Variety of Food and Environmental Samples. 2017, IAFF, Tampa, FL
20. Dutta et al. Performance Evaluation of a Real-time PCR for the Detection of *Cronobacter* spp. in Powdered Infant Formula. 2018. IAFF, Salt Lake City, UT.
21. Dutta et al. Performance Evaluation of a Real-time PCR for the Simultaneous Detection of *Salmonella* and STECs in Co-enriched and Wet Pooled Green Leafy Produce. 2018. IAFF, Salt Lake City, UT.

*Technical review provided for scientific accuracy and does not indicate recommendation or endorsement. The opinions and conclusions presented in this article are the sole position of the authors.

Article originally published on November 12, 2018 by Food Safety Tech Mag and can be found at:

https://foodsafetytech.com/feature_article/pcr-or-lamp-food-safety-considerations-when-choosing-molecular-detection-methods/